


# Choosing the Right Gel Imaging System




## APPLICATIONS

- DNA/RNA gels
- Visible protein gels



## USERS

- User friendly software
- Reproducible results




## VIEWING AREA

- Consider your largest gel size
- Zoom lens




## IMAGE QUALITY

- Uniform illumination
- low background



## CAMERA

- Camera resolution
- Pixel size



## FOOTPRINT

- Laboratory space
- All in one system



## COST

- Meets core requirements
- Reliable performance



## CONNECTIVITY

- Storage options (USB or Network)
- Export formats (TIFF and JPEG)